

# PREPENGINE MINI

The PrepEngine Mini was developed through collaborations with leaders in the pharmaceutical and biotech industry to simplify and speed up the content uniformity sample preparation process. The result - a compact benchtop design that incorporates 2 individually powered, high-output stirring stations which can prepare content uniformity samples up to 90% faster than conventional methods.

- **ULTRA-FAST SAMPLE PREPARATION**

Ultra-fast preparation of biotech, cannabis and pharmaceutical samples including; assay, impurity and content uniformity applications. Up to 90% faster than conventional methods.

- **2 INDIVIDUAL STIRRING STATIONS**

Individually activated, high-power stirring stations maximize productivity and significantly increases sample prep throughput.

- **ADJUSTABLE SPEED & TIME**

The speed is adjustable from 500 to 6000 rpm and the time from 5 seconds to 4 hours.

- **RELIABLE & CONSISTENT RESULTS**

From test to test, the sample results are fast, reliable and consistent.

- **PREPTUBE FOR SINGLE OR MULTI-USE**

Preptubes are a clear polypropylene with SS blade and include a polypropylene septa & cap. Available sizes include 50 ml, 250 mL and 500 mL. Also available in Amber for light sensitive dosage forms.

To eliminate the possibility of carry over and cross contamination, Distek recommends single-use; however, a single PrepTube for wet extraction can be used a maximum of five (5) x 10 minute tests.\*

\*All tests are application specific.

- **CANNABIS SAMPLE PREPARATION**

Dedicated sample preparation station for dry homogenizing, dry grinding & wet extraction of cannabis, hemp and other plant material.

Low-cost, single-use PrepTubes designed specifically for dry homogenizing & grinding.



# PREPENGINE MINI

## SPECIFICATIONS

<b>Stirring Stations</b>	2 Individually Powered, High-Output Stirring Stations with LED Status Indicator	
<b>Adjustable Speed</b>	500 to 6000 RPM	
<b>Free Running Speed</b>	500 RPM $\pm$ 100 RPM / 3000 RPM $\pm$ 150 RPM / 5000 RPM $\pm$ 200 RPM	
<b>Stirring Direction</b>	Clockwise / Counterclockwise Stirring Capabilities if the Clockwise Motion has Stalled	
<b>Adjustable Time</b>	5 Seconds to 4 Hours	
<b>Available PrepTubes</b>	Tube Size (mL): Tube Color: Blade Material: Sterility: Cannabis:	50, 250 and 500 Clear (Standard) or Amber Stainless Steel (Standard) or Polypropylene Non-Sterilized (Standard) or Sterilized Clear, Single-Use for Dry Homogenizing & Grinding
<b>PrepTube Material</b>	Polypropylene	
<b>Display</b>	Backlit LCD (Monochrome)	
<b>Dimensions</b>	10" (w) x 5" (h) x 12" (d) / (25 cm x 13 cm x 30 cm)	
<b>Weight</b>	12 lbs. / (5 kg)	
<b>Electrical Power</b>	100-240VAC / 50-60 Hz / 13A (Operating Voltage Pre-Set at Factory)	



250 mL and 500 mL  
(Clear) PrepTubes



Control Panel (Close Up)



888.2.DISTEK  
[distekinc.com](http://distekinc.com) • [info@distekinc.com](mailto:info@distekinc.com)  
 121 North Center Drive • North Brunswick, NJ 08902