

Product sheet

•	Versatile, precise and reliable	•	•	•	•	•	•	•
•	Total Nitrogen analysis with an						•	
•	extensive application range	•	•	•	•	•	•	•



The TSHR Total Nitrogen Analyzer, model TN 6000, is able to detect precise and accurate trace level nitrogen concentrations in a wide range of liquids, solids and LPG/Gas sample types. The analyzer is designed to run routine total nitrogen analysis down from low ppb to high ppm levels.



The sample is introduced by a fully integrated automatic boat/syringe driver, into a heated oxygen free environment to ensure a complete vaporization of the sample. The carrier gas ensures that the vaporized sample will transfer into the combustion zone where oxygen will be added to complete the oxidation of the sample. After the dual zone combustion stage, the combustion gasses flow through a perma-pure dryer tube where all moisture and other potential interferences are removed.

The conditioned combustion gasses will flow towards the chemiluminescence detector where supplied ozone will react with NO to form an excitation stage of NO₂*. The emitted light during the decay will be detected by the photomultiplier tube. The light intensity is directly proportional to the total nitrogen concentration present in the sample.

The TN 6000 can be optional extended with a TS-UVF module and/or TX module to have a complete versatile combustion elemental analysis solution.

Key advantages

Robust and modular design

Precise, Fast and Reliable Total Nitrogen Data

Optional HR 7000 Liquid autosampler for high sample throughput

Boat cooling option for challenging sample matrices



Analytical specifications

TN Liquids TN Solids

Sample introduction Syringe Liquid module Boat Solids module

Working range 0,03 – 10.000 mg/kg 0,5 – 10.000 mg/kg

Sample matrix Light hydrocarbons Heavy hydrocarbons, solids

 Quantity of Sample
 1 – 100 uL
 0,1 – 100 mg

 Analysis time
 3 - 6 minutes
 4 – 10 minutes

 Relative Standard Deviation*
 < 2% (> 1 ppm)
 < 5% (> 1 ppm)

Regulatory Compliance ASTM D4629, ASTM D5762, ASTM D6069, ASTM D7184, UOP 971, UOP 936

*Depend on typical application and sample matrix

Technical specifications

Furnace Voltage 2 x 24 V , 50/60 Hz

Furnace Power $2 \times 300 \, \text{W}$ Furnace Temperature Sensor $2 \times \text{Ni-Cr/Ni}$

Furnace configuration Dual temperature controlled

Furnace Temperature 1250 °C Max

Type of Analysis Total Nitrogen (TN)

 Detection Principle
 Chemiluminescence (CLD)

 Dimensions
 960 x 390 x 590 mm (WxHxD)

PC operating system Windows 7 or higher

Computer Intel Core i3 / AMD Phenom or better

Software Athena

Optional Supply HR 7000 Liquid Autosampler for 2 mL vials, GM 7000 LPG / Gas Module

Total Sulfur (UV-Fluorescence) detection, Total Chlorine (Microcoulometric) detection

Facility requirements

Voltage 115/230 V , 50/60 Hz

Power 1200 W

Gas connector 1/8" swagelok

Gasses O₂ (99,6%) medical grade 2.6 or

O₂ (99,995%) 4.5

Ar (99,998%) technical grade 4.8


Gas pressure 2 − 3 Bar (30-45 psi) Ambient temperature 5 − 35 °C (41 − 95 °F)

Contact info

For Technical support: technical support@tshrinternational.com

For Sales: sales@tshrinternational.com
For all other queries: info@tshrinternational.com
TSHR website: www.tshrinternational.com

For more information contact your local dealer:



your partner

in combustion

elemental analysis



