Fast, reliable and precise Total

Sulfur analysis with enhanced

performance down to ppb levels

The TSHR Total Sulfur Analyzer, model TS 7000, is the ideal solution for both demanding and routine applications with excellent uptimes and delivers accurate total sulfur data. The TS 7000 incorporates a high-end pulsed UV-Fluorescence detector which provides superior stability, low detection limits and exceptional linearity in compliance with stringent regulatory methods.

The sample is introduced by syringe Liquid autosampler, model HR 7000, in a heated oxygen free environment to ensure a complete vaporization of the sample. The carrier gas ensures that the vaporized sample will carry into the combustion zone where oxygen will be added to complete the oxidation of the sample. After the dual zone combustion stage, the gasses go through a permapure dryer tube where all moisture and other potential interferences are removed.

The conditioned combustion gasses will flow towards the UV-Fluorescence detector where a Xenon flashlamp pulsates UV light which excites the ${\rm SO_2}$ molecules. During decay to lower energy states, UV light is emitted, and detected by the photomultiplier tube. The UV light emission intensity is directly proportional to the concentration of total sulfur present in the sample.

The HR 7000 liquid autosampler is designed for automatic sample throughput at high performance and reliability. This solution provides a most compact fully automated total sulfur analyzer, which can be extended to a capacity of 121 samples automatic analysis or having heated sample tray capabilities.

The TS 7000 analyzes low and high boiling liquids as well as LPG/Gas samples and fully complies with ASTM, DIN, IP and CEN regulatory methods.

Key advantages

Compact and Robust design

Fast, Precise and Reliable Total Sulfur Data

Fully Automated Analysis by HR 7000 Liquid Autosampler

High performance pulsed UV-Fluorescence detector

Enhanced application range for liquids and LPG/Gasses



Analytical specifications

Sample Matrix* Liquid Organics

Working range 0,03 – 10000 mg/kg

Quantity of Sample 1-80 uLAnalysis time 3-6 minutesRelative Standard Deviation* <2% (> 1 ppm)

Type of sample

High & Low boiling point sample

Highest boiling point

450 deg C (subject to sample matrix)

Regulatory Compliance ASTM D5453, ASTM D6667 ASTM D7183, ASTM D7551, ISO 20846, UOP 987

*Depend on typical application and sample matrix

Technical specifications

Furnace Voltage $2 \times 24 \text{ V}$, 50/60 Hz

Furnace Powe $2 \times 300 \, \mathrm{W}$ Furnace Temperature Sensor $2 \times \mathrm{Ni\text{-}Cr/Ni}$

Furnace configuration Dual temperature controlled

Furnace Temperature 1250 °C Max

Type of Analysis Total Sulfur (TS) (optional Total Nitrogen chemiluminescence)

Detection Principle UV-Fluorescence

Dimensions 600 x 1100 x 580 mm (WxHxD)

PC operating system Windows 7 or higher

Computer Intel Core i3 / AMD Phenom or better

Software Athena

Standard Supply** HR 7000 Liquid Autosampler for 2 mL vials

Optional Supply GM 7000 LPG / Gas Module

** HR 7000 model Liquid Autosampler need to be selected for operation of TS 7000

Facility requirements

Voltage 115/230 V , 50/60 Hz

Power 1200 W

Gas connector 1/8" swagelok

Gasses O_2 (99,6%) medical grade 2.6

or O₂ (99,995%) 4.5

Ar (99,998%) technical grade 4.8

Gas pressure 2-3 Bar (30-45 psi) Ambient temperature 5-35 °C (41 -95 °F)

Contact info

For Technical support: technical support@tshrinternational.com

For Sales: sales@tshrinternational.com
For all other queries: info@tshrinternational.com
TSHR website: www.tshrinternational.com

For more information contact your local dealer:



your partner

in combustion

elemental analysis



