
Reliable, fast and accurate

Total Sulfur analysis with enhanced

application versatility

The TSHR Total Sulfur Analyzer, model TS 6000, is the ideal solution for both research and routine applications with excellent modularity and delivers accurate total sulfur data. The TS 6000 incorporates a high-end pulsed UV-Fluorescence detector which provides superior stability, low detection limits and exceptional linearity in compliance with stringent regulatory methods.



The sample is introduced by a fully integrated automatic boat/syringe driver, into a heated oxygen free environment to ensure a complete vaporization of the sample. The carrier gas ensures that the vaporized sample will transfer into the combustion zone where oxygen will be added to complete the oxidation of the sample. After the dual zone combustion stage, the gasses go through a permapure dryer tube where all moisture and other potential interferences are removed.

Key advantages

High performance pulsed UV-Fluorescence detector

Robust and modular design

Boat cooling option for challenging sample matrices

Optional HR 7000 Liquid autosampler for high sample throughput

The conditioned combustion gasses will flow towards the UV-Fluorescence detector where a Xenon flashlamp pulsates UV light which excites the SO_2 molecules. During decay to lower energy states, UV light is emitted, and detected by the photomultiplier tube. The UV light emission intensity is directly proportional to the concentration of total sulfur present in the sample.

The TS 6000 can analyze low and high boiling liquids, solids and LPG/Gas samples and fully complies with ASTM, DIN, ISO, IP and CEN regulatory methods.



Analytical specifications

TS Liquids TS Solids

Sample introduction Syringe Liquid module Boat Solids module

Working range 0,05 – 10.000 mg/kg 0,5 – 10.000 mg/kg

Sample matrix Light hydrocarbons Heavy hydrocarbons, solids

Quantity of Sample1-100 uL0,1-100 mgAnalysis time3-6 minutes4-10 minutesRelative Standard Deviation*< 2% (> 1 ppm)< 5% (> 1 ppm)

Regulatory Compliance ASTM D5453, ASTM D6667, ASTM D7183

*Depend on typical application and sample matrix

Technical specifications

Furnace Voltage 2 x 24 V , 50/60 Hz

Furnace Power $2 \times 300 \, \text{W}$ Furnace Temperature Sensor $2 \times \text{Ni-Cr/Ni}$

Furnace configuration Dual temperature controlled

Furnace Temperature 1250 °C Max

Type of Analysis Total Sulfur (TS)

Detection Principle UV-Fluorescence

Dimensions 960 x 390 x 590 mm (WxHxD)

PC operating system Windows 7 or higher

Computer Intel Core i3 / AMD Phenom or better

Software Athena

Optional Supply HR 7000 Liquid Autosampler for 2 mL vials, GM 7000 LPG / Gas Module

Total Nitrogen (Chemiluminescence) detection, Total Chlorine (Microcoulometric) detection

Facility requirements

Voltage 115/230 V , 50/60 Hz

Power 1200 W
Gas connector 1/8" swagelok

Gasses O₂ (99,6%) medical grade 2.6 or

O₂ (99,995%) 4.5

Ar (99,998% technical grade 4.8


Gas pressure 2-3 Bar (30-45 psi) Ambient temperature 5-35 °C (41 -95 °F)

Contact info

For Technical support: technical support@tshrinternational.com

For Sales: sales@tshrinternational.com
For all other queries: info@tshrinternational.com
TSHR website: www.tshrinternational.com

For more information contact your local dealer:



your partner

in combustion

elemental analysis

