BG Series

Belt Grinders

The BG-Series Wet/Dry Belt Grinders from LECO are available in both single and dual belt models and are engineered for durability, safety, and convenience.





Electric Drive

Heavy-duty 1 hp (0.75 kW) motor provides you with superior spectrographic and metallographic sample preparation.

Improved Belt-Tension and Cleaning

The BG30 and BG32 both feature an enclosed belt tension mechanism that effectively eliminates belt slippage and improves drain flow. Debris removal is facilitated by means of a convenient clean-out access area, and clogging is minimized as a result of improved water flow to the belt.



Safety guards on either side of the grinding area to maximize operator safety.

Durability

The BG Series is built to last with high-quality components and corrosion-resistant aluminum on all external surfaces.

Wet or Dry Grinding

The LECO BG30 and BG32 can be purchased individually or combined with the redesigned front-load support stand and either the 20-gallon stainless steel recirculating system or the dust collector/stand.



BG30 on support stand



BG32 on dust collector

