

MSX250 A-Series

Sectioning Machine

Instrument Profile

These versatile sectioning machines offer both linear and radial feed configurations in efficient benchtop models.

Supports up to 12 in (300 mm) blade @ 5.4 hp (4 kW) with integrated electronic braking system

Sealed spindle design for optimal speed, precise sections, and trouble-free operation

Easily accessible sealed electronic enclosure

LED work light

Radial arm with manual locking feature for table sectioning

Channeled coolant return for high flow rates and minimal debris build-up

MSX250A1

Convenient easy-to-operate control panel

Capable of storing up to 10 applications

Automatic pulse feed rate with load control maximizes sectioning efficiency without overheating samples

13.8 in (350 mm) x 9.8 in (250 mm) T-slotted table for flexible sample fixturing; table vise included

LECO
EMPOWERING RESULTS