


# Miscellaneous Accessories


## Writing Boards & Tackboards



**Black Chalkboard**  
**F-5336-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5337-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5338-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---


Surface: washable Duracite  
 Frame Material: anodized aluminum  
 Trough Material: anodized aluminum



**White Magnetic Markerboard**  
**F-5342-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5343-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5344-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---


Surface: 1/4" white magnetic  
 Frame Material: grey mylar vinyl



**Green Chalkboard**  
**F-5333-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5340-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5341-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---

Surface: washable Duracite  
 Frame Material: anodized aluminum  
 Trough Material: anodized aluminum




**Cork Tackboard**  
**F-5418-24** 18"H x 24"L  
**F-5424-36** 24"H x 36"L  
**F-5436-48** 36"H x 48"L

---

Surface: Resilient natural cork  
 Frame Material: anodized aluminum

## Kemresin Balance Table




30" High x 24" Deep x 35" Long  
**K7-302435-BK** Black  
**K7-302435-GR** Grey  
**K7-302435-PT** Putty  
**K7-302435-SL** Slate

---

Work Top: 2" thick Kemresin  
 Work Top Size: 2"T x 24"D x 35"L  
 Pedestal Legs: 2" thick Kemresin  
 Total Weight: 300 pounds

## Marvel Refrigerators & Freezers



**L-110** Refrigerator/Freezer  
 Size: 34<sup>1</sup>/<sub>2</sub>" H x 24" D x 23<sup>7</sup>/<sub>8</sub>" L  
 6.1 cubic feet  
 Door Swing: Right Hand  
 (User changeable with swing reversing kit)  
 Shelves: (3) Adjustable  
 Service: Single-phase,  
 60 hertz, 120 volt AC

---

**L-105** Refrigerator  
 Same as above except  
 without freezer compartment.

---


**L-115** Freezer  
 Same as above except  
 entire unit is a freezer.

## Gillis Stainless Steel Cart



**F-5349-00**  
 Overall Size: 33<sup>3</sup>/<sub>8</sub>"H x 19"D x 31"L  
 Shelf Size: (3) 18"D x 27"L  
 (1" high lip on 3 sides of  
 top & middle shelf)  
 Capacity: 400 pounds  
 Casters: (4) Black rubber  
 4" Diameter – Swivel

## Lund Key Cases



**F-5332-00**  
 Height: 11"  
 Depth: 2<sup>1</sup>/<sub>2</sub>"  
 Length: 17<sup>5</sup>/<sub>16</sub>"  
 Key Capacity: 90 Keys

---

**F-5332-01**  
 Height: 24<sup>3</sup>/<sub>4</sub>"  
 Depth: 2<sup>1</sup>/<sub>2</sub>"  
 Length: 17<sup>5</sup>/<sub>16</sub>"  
 Key Capacity: 200 Keys